

Features

- **Single dose concept –easy of Handling.**
- **Range of profile covered-Thyroid , Fertility, Hepatitis, Torch, Tumor markers, Allergy, Covid Ab .**
- **Has the electronic calibration, factory set**
- **No consumables except the kit (Calibration/Controls included)**
- **No daily/weekly maintenance required**
- **Easy to operate**
- **Maximum M.T.B.F (more than 1100 days)**
- **12 tests/Run -MiniVIDAS -30 tests/Run-VIDAS PC.**