Features

- Discrete, Random Access, Fully Automated, Compact, Bench Top Model
- Dual Probe, Individual Probes for Sample and Reagent with Liquid Level Sensor
- Onboard Cooling for Reagent Tray
- Onboard Laundry
- Minimal Water Consumption, less than 2 Litres per Hour
- Permanent Optical Glass Cuvettes
- Non Vacuumised Waste/Wash Bottles
- Facility of Post-Dilution and Pre-Dilution Run
- Q.C. with L.J. Plot and Westgard Rules
- Dedicated Onboard Mixer
- STAT Sample Capability
- Built in Computer