Features

- Accommodates different microplate geometries U, V & Flat
- Plate shaking facility with variable speed and duration
- Ability to read plate within seconds
- Digital light control with unique circuit for long lamp life
- Bichromatic readings
- Flexible controls and calibrators positioning
- Storage of calibration curves & controls
- External printer connectivity
- Patient report with PID & name
- Various Graph curve fit mode: Linear, Point to Point, 4 parameter log regression, Cubic Spline,
 Polynomial with scales such as Log Abs / Conc, Logit Abs / Log Conc, Abs / Conc, Log Abs / Log Conc etc.
- Matrix mode for result printing
- Storage of plate results
- Data transfer to computer with easy view of results in Excel format